

PIEZO SURGERY

The 3rd Generation

- ◆ New technology
- ◆ Eight optimizations
- ◆ Rhinoplasty function



CHOOSE US

EIGHT OPTIMIZED



New model

vs



Old model

RHINOPLASTY MODE AVAILABLE



◀ OLD MODEL

NEW MODEL ▶

With Rhinoplasty function,
equiped with Rhinoplasty
tips(Optional)



UPGRADED PERISTALTIC PUMP



◀ **OLD MODEL**
water inlet pipe with
1-point force,
less efficiency and
higher noise.



◀ **NEW MODEL**
water inlet pipe with
3-point force,
More efficiency and
lower noise.

OPTIMIZED THE HANDPIECE LIGHT

OLD MODEL ►

the light is still on
after stop working




NEW MODEL ►

the light is off when
stop working



OPTIMIZED THE HANDPIECE



◀ **OLD MODEL**
Internal water inlet pipe

◀ **NEW MODEL**
External water inlet pipe

- Easier to change water inlet pipe
- Disposable water inlet pipe can be used
- No damage to the water inlet pipe during sterilization


IMPROVED THE HANDPIECE HOLDER

OLD MODEL ►



NEW MODEL ►

Handpiece holder is placed on the saline rod, which is more convenient and save space.



OPTIMIZED WORKING MODE



◀ **OLD MODEL**
Need to Choose Tips



◀ **NEW MODEL**
Fully automatically,
no need to choose tips

ADDED COOLING VENTS


OLD MODEL ▶

No cooling vents



NEW MODEL ▶

Added cooling vents
with waterproof design



2 IN 1 FOOT PEDAL



◀ OLD MODEL

Wireless foot pedal



◀ NEW MODEL

Two options for
foot pedal

■ Option 1: wireless foot pedal

■ Option 2: cable foot pedal



www.baru-dent.com