



IMPLANT WISE PRO

◆ Using precise torque calibration to ensure treatment safety.

AI Artificial Calibration



Real-time torque display



◆ Uses 360° visibility HD non-reflective LCD screen



GERMAN BEARING



- ◆ German bearing ensures extremely stable torque output and silent treatment experience.

SWISS MOTOR

- ◆ Swiss motor offers 100-40000 rpm speed range and 5-80 N/cm torque, which can be used in all kinds of treatment and make cutting bone more efficient and accurate.



- ◆ Automatic LED lighting, brightness can be adjustable, and process visualization

CONTRA ANGLE



- ◆ Unique locking structure, easily disassemble, easy to sanitize the interior.
- ◆ German bearing, radial runout $\leq 0.02\text{mm}$, can withstand 80N·cm large torque output.

GUM RINSER



- ◆ Running silently, offering steady water flow.
- ◆ Easily operate and setting, compact and seamless closure, elegant exterior design.

PERISTALTIC PUMP

- ◆ Water inlet pipe with 3-point force, More efficiency and lower noise.



COOLING VENT DESIGN

- ◆ Added cooling vents with waterproof design.



2 IN 1 FOOT PEDAL

- ◆ Wireless foot pedal: signals received within 10 meters.
- ◆ 4 Hours charging time, and standby more than 50 hours.
- ◆ Water flow control, Program Changes, Motor direction control, these functions can be controlled completely by the multi-function foot pedal, the aseptic operation is truly realized, and the hands are completely free.
- ◆ Two options for foot pedal



Bluetooth connection



Wire connection

