




Dr. Talal' Endo Kit


Instruments

Endodontic Ultrasonic Tips




T1/TD1

Calcification removal in root canal, removal of cement, exposing and removing broken files.



T2/TD2

Active irrigation and irrigation activation.




T3/TD3

Access cavity opening and refinement, canal orifice detection, removing calcifications, loosening crowns and posts.




T4/TD4

Pulp stone and base material removal, pulp chamber calcification removal and cleaning.



T5/TD5

Fiber posts removal (coronally), coronal third canal calcification removal.



T6/TD6

Broken file removal, fiber post removal (apically)

